
IKA

designed for scientists

C 7010 / C 7012

Betriebsanleitung	DE	3
Operating instructions	EN	14
Mode d'emploi	FR	25
Advertencias de seguridad	ES	38
Veiligheidsinstructies	NL	40
Norme di sicurezza	IT	42
Säkerhetsanvisningar	SV	44
Sikkerhedshenvisninger	DA	46
Sikkerhetsanvisninger	NO	48
Turvallisuusohjeet	FI	50
Instruções de serviço	PT	54
Wskazówki bezpieczeństwa	PL	56
Bezpečnostní pokyny	CS	58
Biztonsági utasítások	HU	60
Varnostna opozorila	SL	62
Bezpečnostné pokyny	SK	64
Ohutusjuhised	ET	66
Drošības norādes	LV	68
Nurodymai dėl saugumo	LT	70
Указания за безопасност	BG	72
Instrucțiuni de securitate	RO	74

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DE

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DECLARATION OF CONFORMITY

EN

We declare under our sole responsibility that the product to which this declaration relates is in conformity with the Pressure Vessel Directive 2014/68/EU and the regulations AD 2000.

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In order to be able to use the appliance properly and safely, every user must first read the operating instructions and observe the safety instructions contained therein. Take care of these operating instructions and keep them in a place where they can be accessed by everyone.

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Warranty

You have purchased an original IKA laboratory machine which meets the highest engineering and quality standards.

In accordance with IKA guarantee conditions, the guarantee period is 12 months. For claims under the guarantee please contact your local dealer. You may also send the machine direct to our works, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

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Explanation of symbols

This symbol identifies information that is **absolutely essential for your safety and health**. Failure to observe this information may result in negative effects for your health and injuries.

This symbol identifies information that is important for problem-free operation as well as for working with the decomposition vessel. Failure to observe these instructions may result in inaccurate measurement results or damage to the decomposition vessel.

Contents

EN

Warranty	2
Explanation of symbols	14
Contents	15
Safety instructions	16
Device description	18
Unpacking	18
Included with delivery of C 7010 / C 7012	18
Coding of the decomposition vessel	18
Opening the decomposition vessel	19
Preparing the decomposition vessel	19
Closing the decomposition vessel	19
Combustion procedure	20
Cleaning the decomposition vessel	20
Halogens and decomposition	20
Care and maintenance	21
Leakage test	22
Degase	22
Decomposition vessel identification	22
Accessories	23
Technical data	24
Spare parts list	36
Spare parts diagram	37

Safety instructions

Intended purpose:

The C 7010 and C 7012 decomposition vessels must only be used for determining the gross calorific value of solid and liquid combustible materials in accordance with DIN 51900, BS 1016 T5, ISO 1928, ASTM 5468, ASTM 5865, and ASTM 4809. Operation of decomposition vessels is only permitted in combination with the IKA calorimeter C 7000.

Operating requirements:

The maximum amount of energy input into the decomposition vessel must not exceed **40000 J**. (Select the weight of the sample accordingly). The permissible **operating pressure of 230 bar (23MPa)** must not be exceeded. The maximum permissible **operating temperature** must not exceed **50 ° C**. Do not fill the decomposition vessel too full of the sample. Only fill the decomposition vessel with oxygen up to a maximum pressure of **40bar (4MPa)**. Monitor the adjusted pressure on the pressure reducer.

Perform a check before every combustion to ensure there are no leaks (*see the sections entitled Leakage test*).

Explosive substances:

Many substances tend to combust in an explosive manner (for example because of the formation of peroxide). This may cause the decomposition vessel to burst.

The C 7010 and C 7012 standard decomposition vessels must not be used for examinations on samples that are capable of exploding. It is absolutely essential to use a special high-pressure decomposition vessel to contain the sample in these cases!

Notes on the sample:

Substances of which the combustion behavior is not known must be examined for their combustion behavior before combustion in the decomposition vessel C 7010 and C 7012 (dan-

ger of explosion).

If you are burning **unknown samples**, leave the room or **keep a safe distance** between you and the calorimeter.

Benzoic acid must only be burned in the form of pellets! Combustible dust and powder must be compressed into pellets before combustion. Oven-dry dust and powder such as wood chips, hay, straw, etc. burn in an explosive manner! They must be moistened first! Readily combustible liquids with a low vapor pressure must not be come in direct contact with the cotton thread (for example tetramethyl dihydrogen disiloxan)!

Combustion residue, auxiliary materials

In addition, toxic residues of combustion are possible in the form of gasses, ash or precipitates on the inner wall of the decomposition vessel, for example.

Observe the accident prevention requirements applicable to the activity and the work station. Wear personal safety equipment.

When handling combustion samples, combustion residues and auxiliary materials, the appropriate safety requirements must be observed. The following are examples of substances that may cause dangers:

- corrosive
- easily flammable
- capable of exploding
- contaminated with bacteria
- toxic

Oxygen

When working with oxygen, observe the appropriate requirements.

Danger warning: As a compressed gas, oxygen promotes combustion, supports combustion intensively and may react violently with combustible substances. **Do not use any oil or grease!**

Using a crucible made of stainless steel

When using crucibles made of stainless steel, their condition should be carefully checked after every experiment.

A reduction in the thickness of the material may cause the crucible to burn and may damage the C 7010 and C 7012 decomposition vessel.

For reasons of safety, crucibles must not be used any more after a maximum of 25 combustion procedures.

Specification of the decomposition vessel:

The decomposition vessel is manufactured in accordance with the regulation for pressure vessels 2014/68/EU. This can be recognized from the **CE symbol** with the identifying number of the testing station named. The decomposition vessel is a pressure device of Category III. The decomposition vessel has been subjected to an EC prototype test. The CE declaration of conformity represents our guarantee to you that this decomposition vessel complies with the pressure device described in the EC prototype test certificate. The decomposition vessel has been subjected to a pressure test at a **test pressure of 330 bar (33MPa)** and a leak test with oxygen at 30 bar (3MPa).

 C7010 and C7012 decomposition vessels are **experiment autoclaves** and must be tested by a **professionally trained person** each time before they are used.

An individual application is understood here to mean a series of experiments that are performed under roughly the same conditions in terms of pressure and temperature. Experiment autoclaves must be operated in special chambers (C7000).

Repeated tests

The decomposition vessel must be subject to repeated tests (internal tests and pressure tests) by a **person with professional training**. The intervals between tests must be determined by the operator based on experience, operating manner and the material used in the decomposition vessel.

The declaration of conformity loses its validity if mechanical modifications are made to the experiment autoclaves or if stability can no longer be guaranteed as a result of heavy corrosion (for example holes eaten in it by halogens).

The **threading** on the body of the decomposition vessel and cap screw in particular are subject to a high level of mechanical stress and must therefore be monitored regularly for **wear and tear**.

The condition of the seals must be checked for functionality must be ensured by means of a test for leaks.

(see the sections entitled "Leak test")

Pressure tests and service tasks on the decomposition vessel must only be performed by **persons with professional training**.

We recommend that the decomposition vessel be sent into our factory for inspection and repairs if necessary after either 1000 experiments or after one year or, depending on the application, even sooner than this.

Definition of person with professional training

A person with professional training as defined in these operating instructions is someone

1. whose training, knowledge and experience gained through practical activities ensures that that person will perform the tests in a proper manner.
2. who is sufficiently reliable
3. who is not subject to any instructions in terms of testing activity
4. who is equipped with suitable testing equipment if necessary
5. who can provide suitable proof demonstrating compliance with the requirements listed in 1.

Operating pressure containers

National regulations and laws for operating pressure containers must be observed!

Anyone who operates a pressure container must maintain it in proper condition, must monitor it and perform necessary maintenance and repair tasks without delay, and must take measures appropriate for the circumstances to ensure safety.

A pressure container must not be operated if it exhibits defects that could endanger those working with it or third parties.

You can obtain a copy of the pressure vessel regulation from Carl Heymann Press or Beuth Press.

Device description

Decomposition vessels C 7010 and C 7012 are required for determination of the gross calorific value in the IKA C 7000 calorimeter. The are equipped as described in DIN 51 900 T1 5.1. In addition, these decomposition vessels have integrated temperature sensors which allow the temperature to be measured directly in the vessel.

These decomposition vessels are pressure containers in which combustion may be carried out with oxygen atmospheres. The pressure contain consist of the three pressure-bearing parts which are correspondingly marked. They include:

- bottom part
- cover
- cap screw

The device also features built-in components for electric ignition and or holding the combustion crucible (in accordance with DIN 51 900 T1 5.6): These include:

- electrode
- dish holder

The electrode also assumes the function of the non-returnable valve. The valve is used to both fill the vessel with oxygen (C48) and to empty it. The insides of the decomposition vessel are protected by employments from suitable material.

For easy handling, the bottom part of the decomposition vessel is provided with a sealing (item 18) on which the device can also be set up on to laboratory tables. The integrated board (item 20) which is used for transmission of the temperature signal and for recognition of the decomposition vessel makes it possible to position the decomposition vessel wherever required.

Unpacking

Please unpack the decomposition vessel carefully and check it for damage. It is important that any transport damage is detected when the vessel is unpacked. If necessary, and inventory of the damage should be made immediately (postal service, railway, hou-lage company).

Included with delivery of C 7010 / C 7012

- 1 x decomposition vessel C 7010 or C 7012
- 1 x o-ring item 4
- 1 x cotton thread C7104 (500pieces)
- 1 x o-ring item 7
- 1 x ignition wire C7010 or C7012
- 1 x double open end spanner
- 1 x quartz crucible C4
- 1 x operating instruction C 7010 / C 7012


Coding of the decomposition vessel

When working with the C 7000FZ calorimeter system, a maximum of 8 decomposition vessels can be used. This is made possible through coding of the decomposition vessel from No. 0 to No. 7.

The calorimeter automatically recognizes the decomposition vessel that is presently being used for a measurement.

Coding is carried out by separating the conductor line on the printed circuit board (item 20) located in the bottom part of the decomposition vessel. the conductorlines are marked with I, II and III. (see Figure 1)

It is utmost importance that the thick conductor lines are not damaged when the thin conductor lines are separated. It is advisable to use a sharp knife in order to cut the thin conductor lines at a distance of 2 mm. The piece in between is then removed.

Conductor line I III II

Coding no.:	Conductor line to separate		
	III	II	I
0			
1	X		
2		X	
3	X	X	
4			X
5	X		X
6		X	X
7	X	X	X

Table 1 Coding

Example to illustrate table 1, Coding:

- No.: 0 No conductor line is separated.
 No.: 3 Conductor lines III and II are separated.
 No.: 6 Conductor lines II and I are separated.


Opening the decomposition vessel

- Take hold of the decomposition vessel and firmly grip the bottom part
- Unscrew the cap screw in a counter-clockwise direction to the bottom part
- Use the venting handle in order to pull out the cover. The cover may also be removed by pressing on the sides of the filler neck until the cover is easily pulled out.

The venting handle is part of delivery quantity for the C7000 calorimeter.

Preparing the decomposition vessel

The decomposition vessel is equipped with a fixed ignition wire. A cotton thread has to be fastened to the ignition wire in order to carry out combustion. Fasten cotton thread as shown in the diagram.


Inserting the crucible


The crucible can be placed into the crucible holder after the combustion samples have been weighed.

You can find crucibles among the accessories for the C7010 and C7012

Closing the decomposition vessel

- Put the cover onto the bottom part and press the cover in the bottom part until stop.

- Place the cap screw on to the bottom part and screw handtight.
- The decomposition vessel can now be filled.

Combustion procedure

For the decomposition vessel charging procedure with fuel samples and the general procedure for performing combustion please refer to the relevant standards and the operating instructions for the IKA C 7000 calorimeter system

When working with the C 7000 system, filling of the decomposition vessel with the filling station C 48.

For the C 7000, venting of the decomposition vessel with the venting station C 7030 or with venting botton C 7010.8 supplied with the device is required.

Cleaning the decomposition vessel

After each combustion the condensate has to be cleaned from all of the parts of the decomposition vessel which come into contact with the combustion products. It suffices to wipe out the parts with an absorbent, non-fibrous cloth. Combustion residues in the crucible, e.g. carbon black or ash are removed in the same manner. The decomposition vessel may not be mechanically cleaned under any circumstances if it is not capable of being cleaned in the

manner described (e.g. due to burn-in or pitting).

In such cases please contact your competent service center or send in the decomposition vessel to the factory for inspection and cleaning.

Halogens and decomposition

Pitting which severely affects operation and thus operator safety may be caused on the inner surfaces of the decomposition vessel by combustion products which contain halogen. If the halogen content of the samples is significantly greater than the trace range, then decomposition vessel C 7012, which is specially designed for such applications, must be used.

The C 7010 decomposition vessel may not be used for such applications because of the safety risk and insufficient analysis results.

Decomposition vessel C 7012 is obtainable in a modified form for quantitatively determining halogens and sulfur. This means that the inner surfaces and the fittings of the decomposition vessel are catalytically activated by means of a patented procedure.

Parts thus treated are recognizable because of their opaque, tarnished and partly spotty surfaces. This appearance of this surface is intentional and should not be confused with traces of wear.

This state is important for the following analysis and may not be treated with brushes or other hard materials.

Care and Maintenance

The work described in the following may only be performed by authorized personnel or IKA.

A new leak test is required after any work performed on the fittings of the decomposition vessel.


Furthermore, reference is made to the regulation for pressure vessels 97/23/EC.

Send in the case of service the equipment back in the original packaging. Storage packing is not sufficient for the back dispatch. Use additionally a suitable transportation packing.

Replacing the ignition wire

The C 7010 and C 7012 decomposition vessels are equipped with a fixed ignition wire. Since the ignition wire is subjected to a high degree of material stress during combustion, it should be replaced after about visible wear.

- Unscrew the upper hexagonal nut item 6(IV)/11(III) on the ignition electrode item 5 or the crucible holder item 12 using a double open ended spanner SW 5,5 which is included in the delivery quantity

Replacing the ignition electrode

The ignition electrode has to be replaced if it has become worn as a combustion in the decomposition vessel.

- Open the decomposition vessel as described under „Opening the decomposition vessel“.
- Take hold of the complete cover unit.
- Remove the ignition wire.
- Unscrew the hexagonal nut item 6/11(I) in a counterclockwise direction using the double open ended spanner SW 5,5 while firmly holding the ignition electrode (item 5). Screw the hexagonal nut completely off of the ignition electrode and put it aside for later assembly.
- Pull the ignition electrode out of the cover. Make certain that the insulating bushing (item 9) are removed with the electrode. Clean both parts if required.

A new o-ring (item 7) should be used when the new electrode is inserted.

Information on how to replace the o-ring is provided further on in this chapter. Make certain that there are no dirt particles in the cover when reassembling the unit. Use a paper towel to clean the sealing surfaces.

Installation is performed in the reverse order.

- Put the rate insulating bushing item 8 and 13 and the bearing bushing item 9 on to ignition electrode.
- Insert the new ignition electrode contained in the set of electrodes with deflector (item 2002) into the cover.
- The hexagonal nut (item 6/11(I)) should only be slightly tightened with the spanner.

A leak test should be performed after assembly of the decomposition vessel has been completed.

- Fill the decomposition vessel with oxygen. Use the C48 oxygen station in order to do so.
- Hold the decomposition vessel with the cover upside down in vessel that is completely filled with water. Be careful to only immerse only half of the decomposition vessel at most.

- Check the decomposition vessel for escaping bubbles of oxygen.

The leak test does not replace pressure testing !

Oxygen valve maintenance

The oxygen valve is used to fill and empty the decomposition vessel. It is a non-returnable valve that only closes when the decomposition vessel is full with oxygen.

The ignition electrode item 5 serves at the same time as a valve piston. The sealing surface is in the cover of the decomposition vessel. The O-ring sealing wears through use and because of aggressive combustion gases, which means that the decomposition vessel may not be tight.

Check the sealings for wear in accordance with regular intervals stipulated by the safety provisions. If the decomposition leaks during combustion in the calorimeter the escaping oxygen may damage the decomposition vessel and the calorimeter.

- In order to replace the O-ring the ignition electrode must be removed from the cover as described.
- Remove the O-ring item 7 by first pulling the bearing bushing (item 9) off of the ignition electrode and then removing the O-ring.
- Place the new O-ring from the set of O-rings onto the ignition electrode item 5 from above and then push the bearing bushing item 9 until it again reaches the ignition electrode.

Make certain that the O-ring does not become damaged by sharp edges.

After installation of the ignition electrode the bearing bushing is tightened with the hexagonal nut item 6/11 in order to prevent it from slipping out of place.

Replacing the main sealing

The main sealing item 4 is used to seal off the cover item 2 from the bottom part item 1.

The main sealing and the sealing surface on the bottom part must

be checked for wear by an expert after each combustion operation. The main sealing must be replaced if wear is evident.

- Remove the main sealing (item 4) by using a suitable tool to pry the O-ring out of the groove in the cover (item 2).
- Pull the O-ring toward the side of the cover that is smaller in diameter.

Make certain not to damage the cover with tool. This may cause the decomposition vessel to leak. The sealing that has thus been replaced may not be reused. Replaced the old sealing with a new one from the set of O-rings.

- The new sealing is also mounted on the side of the cover that is smaller in diameter. Place the O-ring on one side the groove in the cover.
- Use your fingers to pull the O-ring over the edge of the cover until it snaps into groove on its own.

Perform a leak test after the decomposition vessel has been assembled.

(see capture "Leakage test")

The leak test does not replace pressure testing !

Replacing the crucible holder

The crucible holder must be replaced if it has become worn as a result of the combustion operations.

- Open the decomposition vessel as described under „Opening the decomposition vessel“.
- Remove ignition wire and ignition electrode as described before.
- Unscrew the hexagonal nut (II) adjacent to the inside cover of the crucible holder in a counterclockwise direction using the SW 5,5 double open ended spanner.
- Screw the crucible holder out of the cover.

Installation of the new crucible holder contained in the set of electrodes with deflector is performed in the reverse order.

The crucible holder has to be aligned before tightening the hexagonal nut. The opening for the crucible must point in the direction of the ignition electrode.

The hexagonal nut is tightened only slightly with the spanner.

Circuit board

The circuit board (item 20) is used for signal transmission. Make certain that the conductor lines are clean.

Leakage test

Perform a leakage test before every combustion process.

Always check the lid is properly sealed with the decomposition vessel after carrying out maintenance.

Follow the safety instructions.

- Fill the decomposition chamber in the calorimeter system or using the C48 oxygen filling station
- Remove the decomposition vessel and immerse fully in a container filled with water
- Check the decomposition vessel for oxygen leaks by watching for bubbles
- If bubbles are seen then the decomposition vessel is not sealed correctly.
- Vent the decomposition vessel and check the seals

Never use a leaking decomposition vessel in the calorimeter.


Degase

Degassing the decomposition vessel takes place manually with the help of the venting handle C 7010.8 in the scope of delivery of the equipment C7000.

The vent knob is supplied for manual release of tension in decomposition vessel.

Escaping decomposition gases might be detrimental to health.


Decomposition vessel identification

The individual parts of the decomposition vessel are clearly identified by means of the engraved serial number and the technical data. The data is in accordance with the regulation for pressure vessels 97/23/EC.

The individual parts of several decomposition vessels (cover, cap screw and bottom part) are not interchangeable.

If one of the parts indicated becomes damaged then the complete decomposition vessel has to be sent to the factory for testing.

Accessories

Ventig handle		C 7010.8
Quartz crucible		C 4
SS-combustion crucible set	(25 pieces)	C 5
Equipment set		C 7050 / C7051
Ignition wire for C7010	(5 pieces)	C 5010.3
Ignition wire for C7012	(2 pieces)	C 5012.3
Cotton threads	(500 pieces)	C 710.4

Gelatin capsules	(100 pieces)	C 9
Acetobutyrate capsules	(100 pieces)	C 10
Benzoic acid (NBS)	(30 g)	C 43
Benzoic acid	(100 g)	C 43A
Benzoic acid tablets	(50 pieces)	C 723
Combustion bags 40x35	(100 pieces)	C 12
Combustion bags 70x40	(100 pieces)	C 12A
Oneway crucible		C 14
Foundation		C 5010.4

Do not use crucibles made of stainless steel for decomposition vessel C 7012

Other accessories can also be found in our **IKA-Catalogue**. The address is included on the back of these operating instructions.


Use of oneway crucible

For the decomposition vessels C 7010 and C 7012 are oneway crucibles C14 ease and shorten the preparation of combustion. They consist of a transparent plastik which burns completely to CO₂ and H₂O.

Advantages:

- The sample can be directly weightet in the oneway crucible.
- The weightof the oneway crucible C14 can be delivered from the weighin-machine to the calorimeter for the calculalte the foreign energy.
- There's no more need for an coton thread.
- The oneway crucible C14 is a burning help. Samples burning incomplete in normal crucibles are burned better with the oneway crucible C14. There by better measuring results are reached.

Replacing the crucible holder


For the use of the oneway crucible C14 the decomposition vessels C 7010 and C 7012 must be put in the in the foundation C5010.4. These foundation C 5010.4 replaces the normal crucible holding.

Take the cover part and until the ignition wire on the crucible holder like discribed before.

Screw down the crucible holder from the shaft and replace it with the foundation C 5010.4.

Installation of the oneway crucible C14

The oneway crucible C14 must be put into the crucible holding of the foundation C 5010.4 so that the ignition wire touches the oneway crucible C14. It would the best to stick it under the ignition wire.

Technical data

Operating pressure PS	max. 230 bar (30MPa)
Test pressure PT	330 bar (33 MPa)
Volume	210 ml
Nuclear temperatures	≥ 1300 °C
Operating pressures	max. 50°C
Sample decomposition	isochor
Oxygen valve	non-returnable valve
	(Vacuum operation not possible)
Crucible holder	as in DIN 51 900
Detection temperature	Pt 100 integratet
Material in contact with the product	

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