2000005297

**IKA®**

## HB 10 basic

## HB 10 digital

## HB 10 control

BETRIEBSANLEITUNG	DE	4
OPERATING INSTRUCTIONS	EN	15
MODE D'EMPLOI	FR	26
使用说明	ZH	38
INDICACIONES DE SEGURIDAD	ES	49
VEILIGHEIDSIINSTRUCTIES	NL	50
NORME DI SICUREZZA	IT	51
SÄKERHETSANVISNINGAR	SV	52
SIKKERHEDSHENVISNINGER	DA	53
SIKKERHETSANVISNINGER	NO	54
TURVALLISUUSOHJEET	FI	55
INSTRUÇÕES DE SERVIÇO	PT	56
WSKAZÓWKI BEZPIECZEŃSTWA	PL	57
BEZPEČNOSTNÍ POKYNY	CS	58
BIZTONSÁGI UTASÍTÁSOK	HU	59
VARNOSTNA OPOZORILA	SL	60
BEZPEČNOSTNÉ POKYNY	SK	61
OHUTUSJUHISED	ET	62
DROŠĪBAS NORĀDES	LV	63
NURODYMAI DĖL SAUGUMO	LT	64
УКАЗАНИЯ ЗА БЕЗОПАСНОСТ	BG	65
INSTRUCȚIUNI DE SECURITATE	RO	66
Υποδείξεις ασφάλειας	EL	67

IKA®-Werke, Germany  
Reg. No. 004343

HB 10 basic


Fig. 1

HB 10 digital/control


Fig. 2

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## Contents

EN

	Page
EU-Declaration of conformity	15
Safety instructions	15
Unpack	16
Correct use	16
Commissioning	17
<b>HB 10 basic</b>	
Operating elements	
Adjustable safety circuit	
<b>HB 10 digital/control</b>	
Heating	
Operating modes	
Safety circuit	
Control of medium temperature	
Interface	
Assembly of spray guard and cover accessories	21
Error message <b>HB 10 basic</b>	22
Error message <b>HB 10 digital/control</b>	23
Maintenance and cleaning	24
Accessories	24
Warranty	24
Technical Data	25

## EU Declaration of conformity

We declare under our sole responsibility that this product corresponds to the directives 2014/35/EU, 2014/30/EU and 2011/65/EU and conforms with the following standards or normative documents: EN 61010-1, EN 61010-2-010, EN 61010-2-051, EN 61326-1, EN 60529 (A1+A2) and EN ISO 12100.

A copy of the complete EU Declaration of Conformity can be requested at [sales@ika.com](mailto:sales@ika.com).

## Safety instructions

### For your protection

- **Read the operating instructions in full before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Wear your personal protective equipment in accordance with the hazard category of the medium to be processed. Otherwise there is a risk of splashing liquids.
- Set up the device in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Prior to each use, always check the device and accessories for damage. Do not use damaged components.
- **CAUTION!** Only process and heat up media that has a flash point higher than the adjusted safe temperature limit of the heating bath that has been set.  
The safe temperature limit of the heating bath must always be set to at least 25 °C lower than the fire point of the media used.
- **DANGER OF BEING BURNT!** When in operation the upper corner of the bath can heat up to 180 °C.
- Prior to filling or emptying the heating bath, the device must be switched off and disconnected from the power supply at the plug.
- Only fill or empty the heating bath when it is cold.
- Empty the heating bath prior to transporting it.
- Never operate the heating bath without tempering medium.
- **CAUTION!** The preferred tempering medium in the heating bath is water (up to approx. 80 °C). Low viscosity silicone oils (50 mPas) with a flashpoint > 260 °C are also permitted. There is a risk of burning when using tempering media with lower flashpoints!
- When using as a water bath, the use of demineralised water is recommended.
- Always observe the minimum fill level of one litre when using oil as the tempering medium.

- Ensure that the IR interface (**HB 10 digital/control**) is not soiled.
- Beware of hazards due to flammable materials.
- Only process media that will not react dangerously to the extra energy produced through processing. This also applies to any extra energy produced in other ways, e.g. through light irradiation.
- Do **not** operate the appliance in explosive atmospheres, with hazardous substances or under water.
- Refer to the operating instructions for the accessories.
- Safe operation is only guaranteed with the accessories described in the "Accessories" chapter.
- Use the HB 10.2 cover or the HB 10.1 spray guard.
- Always disconnect the plug before fitting accessories.
- The appliance **HB 10 basic** does start up again automatically following a cut in the power supply.
- The device is only disconnected from the power supply network if the device power switch is off or the plug is pulled out.
- The socket for the mains cord must be easily accessible.

#### For protection of the equipment

- The voltage stated on the type plate must correspond to the mains voltage.
- Socket must be earthed (protective ground contact).
- Protect the appliance and accessories from bumps and impacts.
- The appliance may only be opened by experts

## Unpacking

### • Unpacking

- Please unpack the device carefully
- In the case of any damage a fact report must be set immediately (post, rail or forwarder)

### • Delivery scope

- <b>HB 10 basic</b>	- <b>HB 10 digital</b>	- <b>HB 10 control</b>
- Heating bath	- Heating bath	- Heating bath
- Mains cable	- Mains cable	- Mains cable
- Manual	- Manual	- Manual

## Correct use

### • Use

- ▶ Tempering the media in the bath container.  
Use with the RV 10 series of rotary evaporators.

### • Range of use

- |                |                |
|----------------|----------------|
| - Laboratories | - Schools      |
| - Pharmacies   | - Universities |

The safety of the user cannot be guaranteed if the appliance is operated with accessories that are not supplied or recommended by the manufacturer or if the appliance is operated improperly contrary to the manufacturer's specifications.

## Commissioning

The unit is ready for service when the mains plug has been plugged in.

Prior to commissioning the device, determine the filling quantity of the tempering medium due to the used evaporator piston size (approx. 2,5 l tempering medium by using a standard 1 l evaporator piston).

### HB 10 basic

Operating elements, see fig. 12

Fig. 12

- **Signal lamp "Power" (G)**
  - ☞ Switch on the device using the main switch (B), see fig. 1.
  - The "Power" signal lamp lights up green.
- **Signal lamp "Heating" (D)**
  - ☞ Use the rotating knob to set the desired temperature of the medium.
  - The "Heating" signal lamp is yellow until the desired temperature is reached.
- **Signal lamp "Error" (C)**
  - ☞ The "Error" signal lamp turns red in the event of an error
  - See the section entitled "Error message"

- **Rotating knob "Temperature setting" (E)**
  - Temperature of medium (area room temperature - 180 °C)
  - ☞ Use the rotating knob to set the desired target value for "Temperature".
  - The "Heating" signal lamp stays lit while energy is being supplied and goes out when the pre-set temperature of the medium has been reached.

### Adjustable safety circuit (F)

The adjustable safety circuit prevents the temperature of the heating bath from getting too high as a result of

- a controller error,
- the "Temperature" rotating knob having been accidentally turned.
- Once the safety temperature has been reached, the device switches off permanently.

Furthermore, dry running of the heating bath is detected.

- The device switches off permanently.

Prior to commissioning, the adjustment knob of the safety circuit "Safe Temp" (F) should be turned all the way to the right (small arrow on the knob pointing to the bottom right).

#### 1. Silicone oils as tempering medium

When increased safety is not necessary, the adjustment knob of safety circuit stays all the way to the right.

Temperature of the tempering medium rises to a maximum of 190 °C in the case of error. Then the safety circuit switches off the heating bath permanently.

#### 2. Demineralised water as tempering medium

When increased safety is not necessary, the adjustment knob of safety circuit stays all the way to the right. As long as there is any demineralised water in the heating bath, the temperature of the tempering medium rises to a maximum of 100 °C in the case of error. If the water is evaporated completely, the temperature of the tempering medium rises to a maximum of 190 °C. Then the safety circuit switches off the heating bath permanently.

### 3. Temperature-sensitive or easily flammable media in the rotary evaporator

When the pre-set bath temperature may not be exceeded under any circumstances set the safety circuit as described below.

#### Setting the safety circuit

- Once the desired temperature of the medium has been reached and the "Heating" signal lamp has blinked three times, the safety circuit can be set.
- Using a suitable tool, slowly turn the adjustment knob (F) counterclockwise until the red signal lamp lights up.
- Then press and slowly turn the adjustment knob (F) clockwise. Once the red signal lamp goes out, turn the adjustment knob (F) a little bit more in a clockwise direction so that the temperature of the safety circuit is somewhat higher than the pre-set target temperature (Turning the adjustment knob (F) by 5° corresponds to an increase in the safety temperature of approximately 8 K).

#### **Reset**

Press the operating knob (F).

## HB 10 digital/control

### Operating elements, see fig. 13


Fig. 13

- ☞ Switch on the device using the main switch (B), see fig. 2.
- ☞ At every start the used software version and the set operating mode are indicated, see fig. 11.

### Setting the target temperature

- ☞ Set the desired target temperature by turning knob (E), LC display, see fig. 8.
- The heating bath temperature is held constant by the control circuit and is additionally monitored by the safety circuit. In the event of an error in the control circuit, the heating bath is switched off permanently by the safety circuit. An error in the control or safety circuit is shown on the display, LC display see fig. 6. The 'Heating' function can no longer be started.
- When switching off the heating function and the medium temperature is higher than 50° C, on the LC display is shown alternately OFF and HOT, see fig. 5 and 7.

## Operating modes

### Operating mode A

- ☞ Switch on the device at the rocker switch (B).
- The heating function is switched off.
- Safe Temp (safety temperature) will appear for 5 secs (fig. 4).
- The target value is set to 20 °C.
- The target value setting can be adjusted.
- The controller tempering medium is set to water.
- The heating function must be restarted following a power failure.
- The following settings are possible/adjustable:
  - the target value,
  - the safety circuit,
  - the controller tempering medium.
- When the device is switched on at the mains switch, A appears (fig. 11).

*Factory setting: Mode A*

### Operating mode B

- ☞ Switch on the device at the rocker switch (B).
- The heating function is switched off.
- Safe Temp (safety temperature) will appear for 5 secs (fig. 4).
- The target value is set to 20 °C or the most recent temperature setting.
- The controller tempering medium is set to water or the most recent medium setting.
- The heating function must be restarted following a power failure.
- The following settings are possible/adjustable:
  - the target value setting,
  - the controller tempering medium.
- The safety circuit is set to the most recently set value and cannot be adjusted.
- When the device is switched on at the mains switch, B appears.

### Operating mode C

- ☞ Switch on the device at the rocker switch (B).
- The heating function is switched on or off depending on the last selected setting.
- Safe Temp (safety temperature) will appear for 5 secs (fig. 4).
- The last target value set in mode B applies.
- The setting controller tempering medium set before the device was last switched off in mode B applies.
- The following settings cannot be adjusted:
  - the target value,
  - the safety circuit,
  - the controller tempering medium.
- When the device is switched on at the mains switch, C appears.

### Switching operating modes

- ☞ Switch the device off at the mains switch.
- ☞ Press and hold knob (E) and switch the device on at the mains switch. Release the knob (E) after approx. 2 seconds.
- Sequence A, B, C, A etc.
- The operating mode setting appears on the display when starting up the device (see fig. 11).

### **Safety circuit**

The adjustable safety circuit prevents the temperature of the heating bath from getting too high as a result of

- a controller error,
- the "Temperature" rotating knob having been accidentally turned.
- Once the safety temperature has been reached, the device switches off permanently.

Furthermore, dry running of the heating bath is detected.

➤ The device switches off permanently.

When the safety circuit is set, an upper temperature limit between 50 and 190° C is defined.

### 1. Silicone oils as tempering media

When increased safety is not necessary the temperature of the safety circuit stays at 190 °C.

Temperature of the tempering medium rises to a maximum of 190 °C in the case of error. Then the safety circuit switches the heating bath permanently off.

### 2. Demineralised water as the tempering medium

When increased safety is not necessary, the adjustment knob of safety circuit stays all the way to the right. As long as there is any demineralised water in the heating bath, the temperature of the tempering medium rises to a maximum of 100 °C in the case of error. If the water is evaporated completely, the temperature of the tempering medium rises to a maximum of 190 °C. Then the safety circuit switches off the heating bath permanently.

### 3. Temperature-sensitive or easily flammable media in the rotary evaporator

When the pre-set bath temperature may not be exceeded under any circumstances set the safety circuit as described below.

#### Adjusting the safety circuit

☞ Set the desired safety temperature by turning knob (E) within 5 seconds of the mask (fig. 4) appearing on the LC display.

#### Checking the safety circuit

The safety circuit has to be checked by the user once a year.

☞ Fill the heating bath with 1 litre of water as tempering medium.

☞ Set the safety temperature to 100° C.

☞ Set the target temperature to 80° C.

☞ Start the 'Heating' function with the ON/OFF key (B).

☞ When the target temperature is reached, turn the heating bath off and on again at the main switch.

☞ Set the safety temperature to 70° C.

➤ The medium temperature is thus 10 K above the safety temperature, the safety circuit trips, the LC display shows E 24, see fig. 6.

### **Control of the medium temperature**

The medium temperature is limited by the set safety temperature. The medium temperature is controlled by means of a PID controller. The medium temperature is measured by the PT 1000 temperature sensor and is heated up as quickly as possible without overshooting to the set temperature.

The PID controller adapts to the different tempering media and ensures optimum temperature control with minimum temperature drift and fluctuation.

☞ Set the desired medium temperature between room temperature and 180° C by turning knob (E). The max. medium temperature setting is the safety temperature minus 10 °C.

➤ LCD display, see fig. 8

☞ Start the 'Heating' function with the ON/OFF key (B).

➤ An animated heating symbol appears on the LC display, see fig. 8.

➤ The heating bath is heated up to the set temperature.

➤ The target and actual temperature referred to the medium are shown on the LC display, see fig. 8.

#### Select tempering medium

➤ If a target value > 90° C is set, the PID controller is adapted specially for oil.

➤ The symbol OIL appears on the LC display, see fig. 9.

➤ If the target value is reset to 20 °C the PID controller is adjusted again specially to water.

➤ The OIL symbol will disappear from the LCD display, see e.g. fig. 8.

## Interface

### Data transfer via IR interface

Fig.14

The heating bath transfers data via IR interfaces. These are located on the left display side of the heating bath or on the right side of the drive unit. Do not place any objects between the two operating units as otherwise the data transfer may be interrupted (fig. 14)!

### Remote Mode

Using the laboratory device software "labworldsoft" and the rotary evaporator RV 10 digital/ control the device can be operated in "Remote" mode. In this mode, the device can no longer be hand-operated.

➤ LCD display see fig. 10

*Note: Please note the system requirements as well as the operating instructions and help section included with the software.*

## Assembly of spray guard and cover accessories

### Assembly spray guard HB 10.1

- ☞ Assemble the spray guard as illustrated.
- Place the bracket of the spray guard on the front edge of the heating bath.
- Attach the spray guard by tightening the screws using a screwdriver.


Fig.15

### Assembly Cover HB 10.2

- ☞ Assemble spray guard and cover as illustrated.
- Push the angled sheet on the back part of the cover into the guide rails at the rear of the heating bath.


Fig.16

- Screw the back and front part of the cover to the spray guard.


Fig.17

## Error message HB 10 basic

Any malfunctions during operation will be identified by lighting up the red signal lamp "Error".

Proceed as follows in such cases:

Error message	Cause	Effect	Correction
<b>"Error" signal lamp lights up</b>	<ul style="list-style-type: none"><li>- Safety circuit temperature set too low</li><li>- Dry run</li><li>- Controller set higher than safety circuit</li></ul>	Device off	<ul style="list-style-type: none"><li>- Let device cool down</li><li>- Reset</li><li>- Check settings and correct if necessary</li><li>- Top up tempering medium</li></ul>
<b>"Error" signal lamp lights up again after reset, also in case of max. adjusted safety temperature limit</b>	<ul style="list-style-type: none"><li>- Faulty controller</li></ul>	Device defect	<ul style="list-style-type: none"><li>- Contact the service department.</li></ul>

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department
- Send the device for repair, including a short description of the fault.

## Error message HB 10 digital/control

Any malfunctions during operation will be identified by an error message on the display (H), see fig. 6.

Proceed as follows in such cases:

- ☞ Switch off device using the main switch (B)
- ☞ Carry out corrective measures
- ☞ Restart device

Error message	Cause	Effect	Correction
<b>E 2</b>	- In remote operation (PC), no communication between RV 10 and HB 10 - IR connection to RV 10 disconnected	Heating off	- Switch off the device - Clean the IR interface or remove objects from around the IR interface - Switch on the device
<b>E 3</b>	- Temperature inside the device is too high	Heating off	- Switch off the device - Let the device cool down - Switch on the device
<b>E 9</b>	- Fault in storing the safety circuit set temperatures- - Memory chip (EPROM) defective	Heating off	- Switch off the device - Let the device cool down - Switch on the device
<b>E 21</b>	- Safety relay does not open	Heating off	- Switch off the device and switch on again
<b>E 23</b>	- Adjustable safety circuit defective	Heating off	- Switch off the device and switch on again - Allow the device to cool down - Switch on the device
<b>E 24</b>	- Set safety temperature exceeded	Heating off	- Switch off the device and switch on again - Allow the device to cool down - Switch on the device again - Check the safety temperature setting
<b>E 25</b>	- The heater control circuit switch (TRIAC) is faulty The heater or the supply line is disconnected	Heating off	- Switch off the device and switch on again
<b>E 26</b>	- Dry running	Heating off	- Let the device cool down - Fill up the tempering medium - Switch off the device and switch on again
<b>E 27</b>	- Error in calibration	Heating off	- Switch off the device and switch on again
<b>E 28</b>	- Breakage of the controller sensor	Heating off	- Switch off the device and switch on again
<b>E 29</b>	- Short-circuit in the safety sensor	Heating off	- Switch off the device and switch on again
<b>E 30</b>	- Short-circuit in the controller sensor	Heating off	- Switch off the device and switch on again

If the actions described fail to resolve the fault or another error code is displayed then take one of the following steps:

- Contact the service department
- Send the device for repair, including a short description of the fault.

## Maintenance and cleaning

The device is maintenance-free.

### Cleaning

For cleaning disconnect the main plug.

Only use cleansing agents which have been recommended by **IKA®**.

To clean the drive use only water with a detergent that contains tensides, or use isopropylalcohol for stubborn soiling.

### Spare parts order

When ordering spare parts, please give:

- Machine type
- Manufacturing number, see type plate
- Item and designation of the spare part, see [www.ika.net](http://www.ika.net).

### Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

## Accessories

- HB 10.1 Spray guard
- HB 10.2 Cover (incl. spray guard)

## Warranty

In accordance with **IKA®** warranty conditions, the warranty period is 24 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

## Technical Data

		HB 10 basic	HB 10 digital/control
Operating voltage range	Vac		220 - 240 ± 10% 100 - 120 ± 10%
Rated voltage	Vac		115, 230 - 50/60 Hz
Frequency	Hz		50 / 60
Power input at nominal voltage	W		1350
Fuses			220-240 V T 10 A (5x20) 100-120 V T 16 A (5x20)
<i>Function Heating</i>			
Heating output at nominal voltage	W		1350
Temperature range medium	°C		ambient temperature - 180 °C
Temperature indication		infinitely	LCD-Display
Heating temperature setting		scale	±1
Setting accuracy heating temperature	K	±5	±1
Controller deviation (3l water/ 90 °C)	K	±5	±1
Controller deviation (3l silicone oil/ 50 mPas/ 150 °C)	K	±5	±2
Adjustable safety circuit	°C	60-190	50-190
Min. volume bath container	l		1
Max. volume bath container	l		3
Materials in contact with medium			stainless steel 1.4301
Safety circuit	°C		adjustable
Dimension (W x D x H)	mm		295 x 265 x190
Weight	kg		2,8
Class designation acc. DIN 12876			II
Protection acc. to DIN EN 60529			IP 21
Per. ambient temperature	°C		5 - 40
Perm. humidity	%		80
Protection class			I
Overvoltage category			II
Contamination level			2
Operation at a terrestrial altitude	m		max. 2000 about sea level

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