



Pre-installation

1 Set up your device

1. **Unpack** your GridION device*
2. Attach the **cables** and **peripherals** (as shown opposite)
3. Connect the **power supply**
4. Press the **power button**

* Place device on a well-supported, strong, clean bench. Allow 30 cm clearance rear and sides, and do not cover ventilation grilles. See User Manual for detailed installation advice.



! Only those ports and connections shaded in blue should be used.

× These ports cannot be used to connect to a monitor.

† If using HDMI only monitor, use the included DisplayPort to HDMI adapter.

2 Log in

1. **Log in** to your GridION Password: **grid**
2. Open **MinKNOW™** Click the wheel icon on the desktop to load MinKNOW, the device operating software
3. **Log in to MinKNOW** Use your Oxford Nanopore account details
4. Complete the **on-screen tutorials** Familiarise yourself with MinKNOW and device operation



3 Perform hardware check

A **hardware check** is required before performing your first GridION sequencing run.

To initiate the hardware check, see the workflow below or follow the on-screen tutorials. You will require your five Configuration Test Cells for this step.




For troubleshooting, refer to the User Manual.

4

Update software

Ensure optimal device and sequencing performance — view and download MinKNOW updates:

 **Host settings**
(in side bar)

 **Software**
(in side bar)

 **Get update**
(in main panel)

Power off the device (step 5), wait 10 seconds, and turn back on to apply the update.


5

Power off

Follow the workflow below to correctly power off your device:

 **Host settings**
(in side bar)

 **Shutdown**
(in main panel)

 **Shutdown**
(in pop-up box)

6

Discover the Nanopore Community



community.nanoporetech.com

Ensure the **success** of your nanopore sequencing project and **stay up-to-date** with the latest technology and **protocol updates**



Additional information



Storage and warranty

GridION device: Do not expose device to temperatures outside of 5–40°C. Flow cell warranty:
community.nanoporetech.com/to/warranty



Recycle used flow cells

Oxford Nanopore is committed to environmental sustainability. Please help us by sending back your MinION™ Flow Cells for recycling. Find out how: community.nanoporetech.com/support/returns



Support

Our technical specialists will be in contact with you shortly to ensure the device is set up correctly. If you encounter issues at any point, please contact support@nanoporetech.com or visit community.nanoporetech.com/support



GridION_{Mk1}

Quick start guide