IKA® C 248

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Contents

	Page
Declaration of Conformity	7
Warranty	7
Warning symbols	7
Safety instructions	8
Correct use and description of the device	8
Unpacking	8
Installation	9
Commissioning	9
Maintenance and cleaning	9
Accessories	10
Technical data	10
List of spare parts-Illustration of spar parts	2

Declaration of conformity

We declare under our sole responsibility that this product corresponds to the regulations 2014/68/EU and conforms with the standards or standardized documents: EN ISO 12100.

Warranty

In accordance with **IKA**[®] warranty conditions, the warranty period is 12 months. For claims under the warranty please contact your local dealer. You may also send the machine direct to our factory, enclosing the delivery invoice and giving reasons for the claim. You will be liable for freight costs.

The warranty does not cover worn out parts, nor does it apply to faults resulting from improper use, insufficient care or maintenance not carried out in accordance with the instructions in this operating manual.

Warning symbols

 DANGER	Indicates an (extremely) hazardous situation, which, if not avoided, will result in death, serious injury.
 WARNING	Indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.
 CAUTION	Indicates a potentially hazardous situation, which, if not avoided, can result in injury.
 NOTICE	Indicates practices which, if not avoided, can result in equipment damage.

Safety instructions

For your protection

Read the operating instructions in full before starting up and follow the safety instructions.

- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.

Work with the device

DANGER

Before use, check the outlet pressure on the pressure reducer and compare it with the pressure displayed on the pressure gauge.

Oxygen as a compressed gas is oxidising; intensively aids combustion; can react violently with flammable materials.

- The filling time should not be set to less than 30 seconds. This will prevent combustible substances from being whirled around in the crucible.
- Insufficient ventilation can cause formation of explosive mixtures. Therefore, only use the device in well-ventilated areas.
- Check the seals for damage before each use.
- After finishing work, the oxygen supply on the pressure reducer must be closed.

For protection of the device

DANGER

Maintenance tasks should be performed without any excess pressure on the system.

- Keep lines and fittings through which oxygen flows free of grease.

Correct use and description of the device

Use

The IKA® C 248 oxygen filling station is only suitable for filling decomposition vessels C 5010 with oxygen.

Description of the device

- The 2 m pressure hose, included with delivery, can easily be laid to connect to the oxygen bottle.
- Therefore, the location of the C 248 filling station does not depend on the oxygen bottle being in the immediate vicinity.
- The filling pressure that is set and the pressure in the decomposition vessel can be read directly from the manometer of the C 248 (11), thus preventing insufficient filling quantities because of not observing the minimum filling period.
- Constant opening and reasoning on the oxygen bottle are void because the automatic non-return valve (bulb compl. (2)) only opens if a decomposition vessel is being filled.

Unpacking

Unpacking

- Please unpack the device carefully.
- In the case of any damage a detailed report must be sent immediately (post, rail or forwarder).

Delivery scope

- C 248 oxygen filling station
- 3 x O-rings
- Screw fitting
- Gasket
- Operating instruction


Installation

Connecting the pressure hose

The pressure hose, included with delivery, is designed for a maximum pressure of 40 bar. It is 2 m long and easily be laid to connect to the oxygen bottle or another oxygen inlet with pressure reducing valve (smallest permissible bending radius 80 mm).

There is an R1/4inch threading (SW 17) on the C 29 reducing valve (accessory) of the oxygen bottle to which the other end of the pressure hose can be connected.

An appropriate adapter is included for US pressure reducers.

Decomposition vessels

The IKA® decomposition vessel C 5010 is provided with an identification on the base or on the cap screw that contains information including the year of manufacturing number.

Commissioning

The required filling pressure -usually 30 bar- can be set on the pressure reducer valve of the oxygen bottle. The adjusted pressure can also be read on the pressure gauge (11).

The bulb compl. (2) integrated into the oxygen filling station prevents oxygen from flowing out when this is not intended.

When the task is complete, the main valve on the oxygen supply should be closed.

The decomposition vessel should be placed on the marked position.

All you need to do then to fill up the decomposition vessel is to swing the rocking arm downwards. The increase in pressure (ca. 2 bar) can be observed on the pressure gauge.

The filling time should not be less than 30 seconds.

Maintenance and cleaning

The C 248 oxygen filling station requires no maintenance to operate. It is subject only to the natural aging of the components and their statistical failure rate.

When ordering spare parts, please indicate the manufacturing number on the rating plate, the device type and the item number and description of the spare parts.

Please send in only units that have been cleaned and that are free of hazardous substances for repairs. Use only cleaning agents approved by IKA® to clean IKA® devices.

To clean the residues below, use the corresponding cleaning agents:

Dye, coloring	Isopropanol
Construction materials	Water containing detergent
Cosmetics	Water containing detergent/isopropanol
Foodstuffs	Water containing detergent Combustible
materials	Water containing detergent

- For materials that are not listed above, please inquire with IKA®.
- Wear protective gloves to clean the units. Electrical devices must not be placed in the cleaning agent to clean them.
- Before using another cleaning or decontamination method than the one recommended by the manufacturer, the user must make certain by consulting with the manufacturer that the intended method will not damage or destroy the unit.

It should be expected that frequent use will result in wear on the O rings in the bulb compl. (2). The filling head (3) is screwed in place and can be unscrewed after the pressure hose is removed. The lower gasket (3007), which is subject to particularly high wear and tear, can be reached without disassembling the pressure hose. Pull with a suitable tool (tweezers) the gasket (3007) down. Put the new gasket into the recess of the decomposition vessel. Place the decomposition vessel under the filling head (3) and fill you the decomposition vessel. The gasket remains at the bulb compl. (2) back.

Repair

In case of repair the device has to be cleaned and free from any materials which may constitute a health hazard.

For this you should use the downloaded print of the „Certificate of compliance“ form from the IKA® website www.ika.com.

If you require servicing, return the appliance in its original packaging. Storage packaging is not sufficient. Please also use suitable transport packaging.

Accessories

- C 29 Pressure gauge, oxygen

Technical data

Permissible operating pressure	bar	max. 40
Pressure gauge	bar	0-50
Typical filling pressure		30
Pressure hose	m	2
Filling speed	sec	> 30
Current filling pressure		readable
Oxygen valve		automatic non-return valve
Dimensions (B x T x H)	mm	130 x 180 x 330
Weight	kg	6

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