## IKA® Calorimeter C 1

### IKA® Calorimeter C 6000 global standards

### IKA® Calorimeter C 6000 isoperibol

Sicherheitshinweise	DE	2	Varnostna opozorila	SL	47
Safety instructions	EN	5	Bezpečnostné pokyny	SK	50
Consignes de sécurité	FR	8	Ohutusjuhised	ET	53
Инструкция по безопасности	RU	11	Drošības norādes	LV	56
Indicaciones de seguridad	ES	14	Nurodymai dėl saugumo	LT	59
Veiligheidsinstructies	NL	17	Инструкции за безопасност	BG	62
Norme di sicurezza	IT	20	Indicații de siguranță	RO	65
Säkerhetsanvisningar	SV	23	Υποδείξεις ασφάλειας	EL	68
Sikkerhedshenvisninger	DA	26	<b>安全说明</b>	ZH	71
Sikkerhetsanvisninger	NO	29	<b>안전 지침</b>	KO	74
Turvallisuusohjeet	FI	32	<b>安全上のご注意</b>	JA	77
Instruções de serviço	PT	35			
Wskazówki bezpieczeństwa	PL	38			
Bezpečnostní pokyny	CS	41			
Biztonsági utasítások	HU	44			

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This "**Safety instructions**" includes only the most important information! It does not replace "**Operating instructions**".

Please read "**Operating instructions**" carefully!

## Explication of warning symbols

**DANGER**

DANGER indicates an imminently hazardous situation, which, if not avoided, will result in death, serious injury.

**WARNING**

WARNING indicates a potentially hazardous situation, which, if not avoided, can result in death, serious injury.

**CAUTION**

CAUTION indicates a potentially hazardous situation, which, if not avoided, can result in injury.

**NOTICE**

NOTICE indicates practices which, if not avoided, can result in equipment damage.

## Safety instructions

- **Read the operating instructions completely before starting up and follow the safety instructions.**
- Keep the operating instructions in a place where it can be accessed by everyone.
- Ensure that only trained staff work with the device.
- Follow the safety instructions, guidelines, occupational health, safety and accident prevention regulations.
- Set up the equipment in a spacious area on an even, stable, clean, non-slip, dry and fireproof surface.
- Do not operate the equipment in explosive atmospheres, with hazardous substances.
- Protect the equipment and accessories from bumping and impacting.
- Check the equipment and accessories for damage before each use. Do not use damaged components.
- Safe operation is only guaranteed with the accessories described in the "**Accessories**" section.
- The socket for the mains cord must be easily accessible.
- Socket must be earthed (protective ground contact).
- The voltage stated on the type plate must correspond to the mains voltage.
- The equipment can only be disconnected from the mains supply by pulling out the mains plug or the connector plug.
- Disconnect the mains plug before attaching or changing any accessories.
- Disconnect the mains plug before cleaning, maintenance and transportation of the unit.
- The equipment must only be opened by trained specialists, even during repair. The equipment must be unplugged from the power supply before opening. Live parts inside the equipment may still be live for some time after unplugging from the power supply.
- Coverings or parts that can be removed from the device without the aid of any tools must be put back on the device again to ensure safe operation, for example to keep foreign objects and liquids, etc. from getting into the device.
- The **IKA®** calorimeter system C1/ C 6000 global standards/isoperibol calorimeter system can be used only for determination of the calorific values of solid and liquid substances to national and international standards (such as **DIN 51900**, **BS 1016 T5**, **ISO 1928**, **ASTM 5468**, **ASTM 5865** and **ASTM 4809**).

**DANGER**

### **C 6000 global standards/isoperibol:**

The maximum extra energy added to the decomposition vessel must not exceed **40,000 J** (select the test mass accordingly). The permitted operating pressure of **230 bar (23 MPa)** must not be exceeded. The maximum permitted operating temperature must not exceed **50 °C**.

**⚠ DANGER** **C 1:** The pressure vessel is integrated into the calorimeter **C 1**. The maximum extra energy added to the pressure vessel must not exceed **40,000 J** (select the test mass accordingly). The permitted operating pressure of **142 bar (14.2 MPa)** must not be exceeded. The maximum permitted operating temperature must not exceed **50 °C**.

**⚠ WARNING** When burning substances containing metals, ensure that the maximum energy input is not exceeded!

**⚠ WARNING** Only fill the decomposition vessel (**C 6010 /C 6012**) or calorimeter system (**C 1**) with oxygen to a maximum pressure of **40 bar (4 MPa)**. Check the pressure setting on the pressure reducer for your oxygen supply. Perform a leakage test before each combustion process.

**⚠ WARNING** Some materials tend to explode when combusted (e.g. due to formation of peroxides), which could cause the decomposition vessel to crack. **The calorimeter can not be used for testing explosive samples.**

**⚠ WARNING** If the burning behavior of a material is unknown, it must be tested before combustion in the decomposition vessel (risk of explosion). If you are burning **unknown samples**, leave the room or **keep enough distance** from the calorimeter.

**⚠ WARNING** Benzoic acid may only be combusted in its pressed form! Flammable dust and powder must be first pressed. Oven-dry dust and powder such as splints, hay, straw etc. explode when combusted! Always wet these materials first! Highly flammable liquids with a low vapor pressure (e.g. tetramethyl dihydrogen disiloxane) must not directly touch the cotton thread! Furthermore, toxic combustion residue in the form of gases, ashes or condensation, for example, is possible in the inner wall of the decomposition vessel.

**⚠ DANGER** Observe the accident prevention requirements applicable to the activity and the work station.

- Wear your personal protective equipment.
- When handling combustion samples, combustion residue and auxiliary materials, please observe the relevant safety regulations. The following materials, for example, could pose a risk:
  - corrosive
  - highly flammable
  - explosive
  - bacteriologically contaminated
  - toxic

**⚠ WARNING** Please observe the relevant regulations when handling **oxygen**. Oxygen as a compressed gas is oxidising; intensively aids combustion; can

**⚠ DANGER** react violently with flammable materials.  
**Do not use oil or grease!**

**⚠ WARNING** **Combustion gases are hazardous to health, therefore the venting hose must be connected to a suitable gas cleaning system or extraction system.**

- Tubes and screwed joints for oxygen must be kept free of grease.
- At the end of the work period, close the main valve for the oxygen supply.
- Perform servicing work only when the equipment is depressurized.

**⚠ WARNING** When using **stainless steel crucibles**, thoroughly check their condition after each experiment. If the material gets thinner, the crucible may catch fire and damage the decomposition vessel. Crucibles must not be used for more than **25** combustions for safety reasons.

**⚠ NOTICE** **C 6000 global standards/isoperibol:** The decomposition vessel is manufactured in accordance with the directive for pressure equipment **97/23/EC**. This is indicated by the **CE symbol** with the ID number of the notified body. The decomposition vessel is a category III pressure device. The decomposition vessel has undergone an EC prototype test. The EC declaration of conformity confirms that this decomposition vessel corresponds to the pressure device described in the EC prototype test certificate. The decomposition vessel has undergone a **pressure test** with test pressure of **330 bar (33 MPa)** and a leakage test with oxygen at **30 bar (3 MPa)**.

**⚠ NOTICE** **C 1:** The EC Declaration of Conformity confirms that this pressure vessel was manufactured by the Pressure Equipment Directive **97/23/EC**. (Article 3, paragraph 3). The pressure vessel has undergone a pressure test with test pressure of **203 bar (20.3 MPa)** and a leakage test with oxygen at **40 bar (4 MPa)**.

**C 6000 global standards/isoperibol:** The CE declaration of conformity confirms that this decomposition vessel corresponds to the pressure device described in the EC prototype test certificate. The decomposition vessel has undergone a **pressure test** with test pressure of **33 MPa** and a leakage test with oxygen at **3 MPa**.

**⚠ NOTICE** **C 6000 global standards/isoperibol:** refer also to the operating instructions for the **C 6010/C6012** decomposition vessel.

Decomposition vessels are experiment autoclaves and must be tested by a **technical expert** after each use.

- Individual use is understood here to include a series of experiments performed under roughly the same conditions in terms of pressure and temperature. Experiment autoclaves must be operated in special chambers

**The Calorimeter system C1 and the decomposition vessels**

(**C 6010/ C 6012**) must undergo repeated tests (internal tests and pressure tests) performed by the **technical expert**. The frequency of these tests is to be determined by the operator on the basis of experience, type of operation and the material used in the decomposition vessel.

- The declaration of conformity becomes invalid if mechanical modifications are carried out to the experiment autoclaves or if tightness can no longer be guaranteed as a result of major corrosion (e.g. pitting by halogens).
- In particular the threads on the pressure vessel **C 1, C 6010, C 6012** are subject to considerable stress and must therefore be checked regularly for wear.

The condition and function of the seals must be checked and ensured by way of a leakage test (**C 6010/C 6012**) or a system (**C 1**) test. **Check the seals for damage before each use.**

If the maintenance, and especially the pressure testing, is not performed or is performed incorrectly, there is a risk the decomposition vessel may burst or an uncontrolled internal fire may occur at the electrodes which could burn away the seals (oxyacetylene torch effect), thus posing a risk to life and limb.

Only **technical experts** can perform pressure tests and service work on the decomposition vessel.

**We recommend that you send the decomposition vessel or the calorimeter system C 1 to our factory for inspection, and if necessary, repair after 1000 tests or after one year or sooner depending on use.**

A corresponding warning message in the display can be turned off by a pressure test and enter an activation code.

**Note:** The device will continue to work!

**C 6000 global standards/isoperibol:** A pressure test must be performed after a decomposition vessel has reached or exceeded the recommended number of ignition cycles. Once the pressure test has been performed a release code can be entered to enable the decomposition vessel to be used for further measurements. The warning message will then disappear!

**Note:** You can then continue to use the device /decomposition vessel.

- Contact the **IKA®** Service Department to perform the pressure test. Comply with the safety instructions in this respect. The user will be given specific instructions on the relevant screen.
- For the purposes of these operating instructions a **technical expert** is someone
  1. who guarantees to conduct the tests properly on the basis of his training, knowledge and experience gained through practical work,
  2. who is sufficiently reliable,
  3. who is not subject to any instructions in terms of the test activity,
  4. who has suitable test equipment if necessary,
  5. who can provide appropriate proof of the requirements listed in 1.
- **National directives and laws must be observed for operating pressure vessels!**
- Anyone operating a pressure vessel must keep it in a proper condition, operate it properly, supervise it, carry out the necessary maintenance and repair work immediately and implement the safety measures required in the circumstances.

Pressure vessels must not be used if they have defects which could pose a risk to staff or third parties.

- You can obtain a copy of the pressure vessel regulation from Beuth Press. (see <http://ec.europa.eu/enterprise/sectors/pressure-and-gas/>).

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