

BLOTTING

Model 785 Vacuum Blotter

- Fast, accurate transfer of DNA or RNA
- Quick, easy setup
- Controlled vacuum pressure for consistent results
- Molded construction for long-term reliability
- Chemically resistant construction

Fast Vacuum Transfer System

The Model 785 vacuum blotter is a nucleic acid transfer apparatus that uses low vacuum to transfer DNA or RNA from agarose gels onto nylon membranes in as little as 90 minutes. This vacuum blotter simplifies every step of the blotting process. Its easy assembly reduces setup times for your experiments, while the vacuum regulator accurately controls vacuum pressure to give you consistent results. The Model 785 vacuum blotter has durable molded construction for long-term, heavy use. Detailed protocols for single-copy gene detection are included in the simple, step-by-step instruction manual. The Model 785 vacuum blotter saves valuable experimental time without compromising results.


The Model 785 vacuum blotter has several innovative features that help to simplify the nucleic acid blotting process and facilitate complete transfer, including:

- Spring latches make assembly and disassembly quick and easy — no screws or nuts are needed to tighten the sealing frame
- Vacuum stage and porous vacuum plate create a uniform vacuum pressure over the entire transfer area, for even transfers and highly reproducible results
- Blank window gaskets are available for transferring differently sized gels
- Easy-to-see green window gasket, combined with simple base unit and sealing frame, make aligning the gel and membrane quick and easy
- Vacuum regulator accurately controls vacuum pressure, whether you use a Bio-Rad vacuum pump or a house vacuum
- Durable, molded construction withstands long-term, heavy use
- Blotter is chemically resistant to commonly used acids and bases (0.25 N HCl or 0.5 N NaOH)
- Tested protocols for overnight-exposure single-copy gene detection are included in the Model 785 vacuum blotter instruction manual




BIO-RAD

Model 785 Vacuum Blotter



Southern transfer. Human genomic DNA (5 µg) was cleaved with (left to right) *Bgl*I, *Hind*III, or *Nhe*I, separated on a 1% Certified™ molecular biology agarose gel, and transferred to a Zeta-Probe® GT membrane using a vacuum of 5" Hg for 90 min. The blot was probed with a ³²P-labeled factor VIII fragment and exposed to X-ray film overnight. Probe provided by the Genetics Institute of Boston.



Northern transfer. Left to right, 1.0, 2.5, 5, 10, or 20 µg of total RNA from CHO HA-1 cells was separated on a glyoxal gel and transferred to a Zeta-Probe® membrane using a vacuum of 3" Hg for 90 min. The blot was probed with a ³²P-labeled hsp70 cDNA fragment and exposed to X-ray film overnight. Probe and RNA provided by Gloria Lee, University of California, San Francisco, CA USA.

Specifications

Dimensions

Overall size (L x W x H)	21 x 12 x 5"
Pump station (L x W x H)	16 x 9 x 8"
Vacuum tubing	1/2" OD x 1/4" ID

Weight

Basic unit	8.5 lb
Vacuum station	9.5 lb
Vacuum regulator	1.5 lb

Ratings

Vacuum pump	20" Hg
Vacuum regulator	0–25" Hg
Gel capacity (maximum)	20 x 25 cm

Required Equipment

Vacuum source	Vacuum pump or house vacuum
Vacuum regulator	Strongly recommended for reproducible results
Liquid trap	Sidearm flask with one-hole stopper (not sold by Bio-Rad Laboratories)
Vacuum tubing	Required for liquid trap assembly

Ordering Information

Catalog #	Description
165-5000	Model 785 Vacuum Blotter With Regulator includes vacuum regulator, base with vacuum stage, porous vacuum plate, reservoir seal O-ring, sealing frame, assorted window gaskets, lid, instructions
165-5001	Model 785 Vacuum Blotter System, 100/120 V , includes vacuum pump, vacuum regulator, base with vacuum stage, porous vacuum plate, reservoir seal O-ring, sealing frame, assorted window gaskets, lid, instructions
165-5002	Model 785 Vacuum Blotter System, 220/240 V
165-5003	Model 785 Vacuum Blotter Basic Unit , includes base with vacuum stage, porous vacuum plate, reservoir seal O-ring, sealing frame, assorted window gaskets, lid, instructions
165-5004	Vacuum Station 100/120 V , includes vacuum pump, vacuum regulator
165-5005	Vacuum Station, 220/240 V
165-5006	Vacuum Regulator
165-5007	Reservoir Seal O-Ring
165-5008	Blank Window Gaskets, 6
165-5012	Porous Vacuum Plate

Ultra Pure DNA Agarose

162-0125	High-Strength Analytical Grade Agarose, 100 g
162-0135	Chromosomal Grade Agarose, 25 g
162-0136	Chromosomal Grade Agarose, 100 g

Blotting Accessories

162-0153	Zeta-Probe Membranes, 9 x 12 cm, 15 sheets
162-0154	Zeta-Probe Membranes, 10 x 15 cm, 15 sheets
162-0155	Zeta-Probe Membranes, 15 x 15 cm, 15 sheets
162-0156	Zeta-Probe Membranes, 15 x 20 cm, 15 sheets
162-0157	Zeta-Probe Membranes, 20 x 20 cm, 15 sheets
162-0158	Zeta-Probe Membranes, 20 x 25 cm, 3 sheets
162-0159	Zeta-Probe Membrane, 30 cm x 3.3 m, 1 roll
162-0165	Zeta-Probe Membrane, 20 cm x 3.3 m, 1 roll
162-0190	Zeta-Probe GT Membranes, 9 x 12 cm, 15 sheets
162-0191	Zeta-Probe GT Membranes, 10 x 15 cm, 15 sheets
162-0192	Zeta-Probe GT Membranes, 15 x 15 cm, 15 sheets
162-0193	Zeta-Probe GT Membranes, 15 x 20 cm, 15 sheets
162-0194	Zeta-Probe GT Membranes, 20 x 20 cm, 15 sheets
162-0195	Zeta-Probe GT Membranes, 20 x 25 cm, 3 sheets
162-0196	Zeta-Probe GT Membrane, 30 cm x 3.3 m, 1 roll
162-0197	Zeta-Probe GT Membrane, 20 cm x 3.3 m, 1 roll
161-0302	SDS (Sodium Dodecyl Sulfate), 1 kg
165-0962	Filter Paper Backing, for stained gels, 35 x 45 cm, 25 sheets
165-0921	Thick Blot Paper, 18 x 34 cm, for Model 224, 443, and 543 slab gel dryers, 25 sheets

BIO-RAD

**Bio-Rad
Laboratories, Inc.**