20000011577

C 200.RC_092016

IKA®

IKA® C 200.RC

Technische Information	DE	4
Ursprungssprache		
Technical information	EN	6
Information technique	FR	8
Información técnica	ES	10
Informacje techniczne	PL	12
Informação técnica	PT	14
Техническая информация	RU	16
技术信息	ZH	18
取扱説明報	JA	20
기술 정보	KO	22

IKA®-Werke, Germany
Reg. No. 004343

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Item	Designation
1	Input nipple, pressure side
2	Output nipple, intake side
3	Quick plug coupling, output
4	Level tube (inner vessel drain)
5	Quick plug coupling, input
6	Hose anchoring (magnet)
7	Plastic coupling
8	T-piece (connection line for the inner vessel drain)
9	Pressure control valve with adjusting cap
10	Connection line, intake side
11	Connection line, pressure side
12	Screw cap for each connection line

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Content

	Page
Device setup	2
Area of application	6
Safety instructions	6
Unpacking	6
Assembly	6
Commissioning	7
Maintenance and cleaning	7
Technical Data	7

EN

Area of application

The **IKA**® C 200.RC filling adapter supplies the C 200 calorimeter automatically and continuously with medium (water). With the RC 2 or KV 600 condenser, you have the option of pumping water at a constant temperature into the tank of the C 200. Manually filling with medium at a constant temperature before the measurement is not required.

Safety instructions

Please read the instruction manual of the condenser and the C 200 calorimeter in full before use and follow the safety instructions.

- Keep the operating instructions in a place where they can be accessed by everyone.
- Ensure that only trained staff work with the appliance.
- Follow the safety instructions, guidelines, occupational health and safety and accident prevention regulations.
- Ensure that the connection lines are connected correctly..
- The condenser used (RC 2 or KV 600) must have an active intake side.
- Remove the device from the mains before assembly/disassembly.
- Before use, check the hose system between the devices for leaks.
- Secure the hoses using the supplied hose anchoring (magnet). Do not bend the hoses
- Check hoses, tubes and bath at regular intervals for possible material fatigue (cracks/leaks).

Unpacking

Unpacking

- Please unpack the device carefully.
- In case of any damage, a detailed report must be sent immediately (post, rail or forwarder).

Delivery scope

- C 200.RC filling adapter (quick plug couplings with level tube)
- Connection line, intake side with T-piece for the inner vessel drain
- Connection line, pressure side with pressure control valve
- Hose anchoring (magnet)

Assembly

- Remove the tank filler on the C 200.

- Press the C 200.RC filling adapter in the tank; the input/output nipple shows the back of the device.


- Insert the connection line, intake side (Fig. 2, (10)) with the quick plug coupling (Fig. 2, (3)) onto the input nipple (Fig. 1, (2)). The coupling engages.
- Screw the other end of the connection line, intake side (Fig. 2, (12)) hand-tight to the condenser connection ((IN), intake side). Shorten the hose, if necessary.

- Insert the connection line, pressure side (Fig. 2, (11)) with the quick plug coupling (Fig. 2, (5)) onto the input nipple (Fig. 1, (1)). The coupling engages.
- Screw the other end of the connection line, pressure side (Fig. 2, (12)) hand-tight to the condenser connection ((OUT), pressure side). Shorten the hose, if necessary.
- Insert the plastic coupling (Fig. 2, (7)) into the C 200 (Standard OUT).

- Secure the hoses in place using the supplied hose holder. The hose holder should be installed 30 mm from the edge of the housing, so that the hoses do not disrupt the area of the cover. To do this, remove the film from the adhesive side and secure the hose holder on top of the C 200. Insert the magnets into the recess and secure the hoses.

Magnet Hose anchoring

Commissioning

Note the operating instructions of the condenser (RC 2 or KV 600) and the C 200 calorimeter.

When the inner vessel is filled with medium during a test, the filling level drops in the condenser. The removed quantity is pumped back into the condenser after the test. Therefore, the medium must not be refilled, otherwise there is a risk that the condenser will overflow.

Procedure:

1. Fill the condenser with medium and set the desired working temperature.
 2. C 200 initial fill; refer to the operating instructions for the device.
 3. Manual control of the flow rate (pressure control valve).
- Fill the condenser up to the max. fill level with a suitable medium.
 - Open the adjusting cap on the pressure control valve: Push the cap upwards and turn it clockwise (+) by 5 turns.
 - Switch on the condenser.
 - Set the working temperature on the condenser (e.g. 22 °C). **IKA® recommends operation at 2800 rpm.**
 - Start the pump function of the condenser so that water is pumped through the pressure hose into the tank of the C 200.
 - Once the condenser has reached the desired temperature, the C 200 switches on.
 - Carry out initial filling with the C 200:
For the initial filling of the outer vessel, the C 200 pumps the medium from the tank of the C 200. To balance the fill level, the medium is pumped from the condenser. This reduces the fill level in the condenser. Refill with water up to the max. fill level. The initial filling takes about 1 minute.
 - The flow rate on the adapter can be adjusted using the pressure control valve.

Reduce the flow rate if the fill level in the fill level indicator of the C 200 is visible and rises. The fill level in the C 200 tank should always be lower than what can be displayed in the fill level indicator of the C 200.

Max. fill level of fill level indicator

1. Turn the adjusting cap on the pressure control valve by one turn counter clockwise (-). Check the fill level, if this falls the correction is sufficient, if it continues to rise, turn it by another turn (-). Repeat the procedure until the medium in the fill level indicator is no longer visible.
2. Press the adjusting cap down to close the pressure control valve.


Maintenance and cleaning

Clean the hoses once per month.

The instrument is maintenance-free. It is only subject to the natural wear and tear of components and their statistical failure rate.

Cleaning

For cleaning disconnect the main plug!

Use only cleaning agents which have been approved by **IKA®** to clean **IKA®** instruments.

Dirt	Cleaning agent
Dyes	Isopropanol
Building materials	Water containing detergent, Isopropanol
Cosmetics	Water containing detergent, Isopropanol
Food	Water containing detergent
Other materials	Water containing detergent

For materials which are not listed, please request information from **IKA®** application support. Wear protective gloves during cleaning the instruments. Electrical instruments may not be placed in the cleansing agent for the purpose of cleaning.

Do not allow moisture to get into the instrument when cleaning. Before using another than the recommended method for cleaning or decontamination, the user must ascertain with **IKA®** that this method does not destroy the instrument.

Repairs:

Please only send devices in for repair that have been cleaned and are free of materials which might present health hazards.

For this, use the "certificate of compliance" form which you can obtain from **IKA®** or can download a version for printing from the **IKA®** website at www.ika.com.

If your appliance requires repair, return it in its original packaging. Storage packaging is not sufficient when sending the device - also use appropriate transport packaging.

Technical Data

Permissible pressure	bar/MPa	2/0.2
Medium		tap water/drinking water quality

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